

# DRAMINSKI 4Vet *mini*

PROFESSIONAL ULTRASOUND SCANNER FOR MIXED PRACTICE VETERINARIANS

- excellent image quality
- operates with different electronic probes
- effective battery
- solid and durable frame body
- possible to connect goggles

also available



WITH GOGGLES



Call:  
+48 89 527 11 30

Write:  
[info@draminski.com](mailto:info@draminski.com)

# DRAMIŃSKI 4Vet mini

PROFESSIONAL ULTRASOUND SCANNER  
FOR MIXED PRACTICE VETERINARIANS



- ◀ It has an excellent image quality which facilitates making a proper diagnosis.
- ◀ It operates with different electronic probes which enlarges a range of tests.
- ◀ It has an effective battery enabling more than 3 hours of continuous work.
- ◀ It has a solid and durable frame body guaranteeing safe work in difficult field conditions.
- ◀ It is possible to connect goggles to the scanner, so that working in direct sunlight becomes more pleasurable.



## DESIGNED TO WORK IN THE FIELD

**Standard package includes:** an ultrasound scanner with a selected electronic probe, 2 external batteries, a charger, micro-USB/USB 2.0 adapter, neck and hip strips, ultrasound gel (250 ml), transport case, user manual.

**Technical specification:**

- Size 23.0 x 17.0 x 6.5 cm,
- Weight 2060 g,
- 7" LED-backlit LCD display.

**Image:**

- B, B+B, B+M
- post-processing – digital image processing
- zoom 60 – 200%,
- freeze,
- regulation: gain, gamma, frequency, focus.

**Probes:**

- convex 2 – 8 MHz, penetration depth up to 30 cm,
- microconvex 4 – 9 MHz, penetration depth up to 15 cm,
- linear 60 mm 5 – 10 MHz, penetration depth up to 12 cm,
- linear 40 mm 6 – 14 MHz, penetration depth up to 6 cm,
- endorectal 4 – 9 MHz, penetration depth up to 12 cm.

**Functions:**

- save images and cine loop with description and date to internal disc,
- transmission of images and cine loop to external data carrier,
- measuring (distance, surface, volume, fetal biometry charts),
- regulation of scanning range.



*Cow – uterus horns*



*Sheep – early pregnancy*